

# Face Visors


## Full Face Protection

The anti-fog, latex free eyeshields are designed to protect your eyes and face by preventing exposure from hazardous or infectious substances.


- Frames available in 10 stylish colours
- Provides overhead, side & front face protection
- One size fits all!
- Optical grade, fog resistant, anti-static replacement shields
- Latex free




### Instruction for Use




1. Insert nosepiece tab into visor slot, so nosepiece slants away from front of visor



2. Mount the facieshield onto the visor pin, then insert the visor's support kuckle through the facieshield slot



3. Repeat step 2 on the opposite side of the visor.



### Order details:

MED1000.05	Full Face visors	pack of 10
MED1000.06	Full Face sheild	pack of 25