

Basket Connection Guide



	SI Digital PC+	SI SA / SI SA 'S'	SI PCF / SI PCF 'S'	PCF PRO	Niagara
General Instruments	Basket & Hose Set: MED1018.1 Hose Set: MP22203	SI SA Basket & Hose Set: MP21853 Hose Set: MP22203 SI SA 'S' Basket & Hose Set: MP21854 Hose Set: MP22203	Basket & Hose Set: MP21832 Hose Set: MP22183 Y Hose Set: MP22199	Basket & Hose Set: MP21852 Hose Set: MP21846	Carrier: MED1437 Baskets only: M20044 / M20043 Basket & Hose Set: MED1425 Hose Set: MP22243 Basket: MP21841 Hose Set: MP22222
Orthopedic Sets			Basket only: MED8603		Carrier: MED1415
SI Endowrist 8mm			Basket & Hose Set: MP22062 Hose FP1: M20066	Basket & Hose Set: MP21907 Hose Set: MP22188	Carrier: MED1455 ¹¹ Or Carrier: MED1426 Hose Set: M20134 Basket & Hose Set: MP21838 ¹ Hose Flush Port 1: MP22212
Si Endowrist 5mm / SS			Basket & Hose Set: MP21833 Hose FP1: M20066 Hose for Shaft Port: MP22201	Basket & Hose Set: MP21907 Hose Set: MP22189	Carrier: MED1426 Hose Set: M20134
Si Stapler			Basket & Hose Set: MP21833 ⁹ Hose FP1: M20066 Hose for Shaft Port: MP18699.10	Basket & Hose Set: MP21907 Hose Set: MP22190	Carrier: MED1426 Hose Set: M20134 Basket & Hose Set: MP21839 ³ Hose Flush Port 1: MP22214 Hose for Shaft Port: MP22215
Xi Endowrist 8mm			Basket & Hose Set: MP21104 ¹⁰ Hose Set: MP21152	Basket & Hose Set: MP22036 Hose Set: MP22018	Carrier: MED1426 Hose Set: M20134 Basket & Hose Set: MED1542 ⁴ Hose Set: MP20966.5 Basket: MP21837 ⁵ Hose Set: MP21133 ⁶
Xi Stapler			Basket & Hose Set: MP21104 ¹⁰ Hose Set: MP21152	Basket & Hose Set: MP21906 Hose Set: MP22018	Carrier: MED1426 Hose Set: M20134 Basket & Hose Set: MP21840 ⁷ Hose Set: MP22217 Basket & Hose Set: MP21836 ⁸ Hose Set: MP21133 ⁶

- 1 Wire basket for 5 Si Endowrists
- 2 Wire basket for 5 Si 8mm Endowrists
- 3 Sheet metal basket for 2 Si Staplers
- 4 Sheet metal basket for 5 Xi Endowrists

- 5 Wire basket for 5 Xi Endowrists
- 6 Only to be used with MP21837 and MP21836 baskets
- 7 Sheet metal basket for 2 Xi Staplers
- 8 Wire basket for 3 Xi Endowrists / 2 Xi Staplers

- 9 Compatible with Si 8mm, Si 5mm, Si Single Site and Si Stapler
- 10 Basket including connector / hose for up to 6 Xi Endowrists or 4 Xi Endowrists and 2 Xi Staplers
- 11 Can only accommodate the Si 8mm Endowrists directly on carrier without the use of baskets